




Grupo Álava

Beyond  
technology

[grupoalava.com](http://grupoalava.com)




# Grupo Álava



## We seamlessly connect manufacturers, system integrators and end users

We identify needs and develop high-tech solutions.





Auscultation &  
Civil Engineering



Communications




Environmental  
technologies




Imaging &  
Photonics



Instrumentation  
& Calibration



Nanotechnology



Predictive



Security



Testing


## Lasers



## Optics and optomechanics



## Spectroscopy solutions: hyperspectral cameras, FTIR, portable spectrometers...



## Systems for optical character.

LASER SAFETY OFFICER




## Laser protection




## Imaging solutions


# Innolas Laser devices: LIBS applications and developments

	<p>▶ OPOs</p>	<p>VIS &amp; NIR emitting OPOs from 1mJ to 150mJ</p>
	<p>▶ Sub-ns Lasers</p>	<p>Active Q-Switch sub-nanosecond pulses with up to 1 J pulse energy</p>
	<p>▶ Diode Pumped Lasers</p>	<p>DPSS Lasers with average power of up to 100 W</p>
	<p>▶ Hybrid Lasers</p>	<p>Diode pumped laser with flash lamp amplification - up to 800mJ at 100 Hz</p>
	<p>▶ Lamp Pumped Lasers</p>	<p>Flash lamp-pumped ns Nd:YAG lasers from 100mJ to 2.5J</p>
	<p>▶ Multi-Pulse Lasers</p>	<p>Two laser sources in one housing</p>



- Analysis of heavy metals contained in soil samples: Changchun downtown Changchun Railway Station (dusty), Nanhu Park (less tourists), Jingyuetan Park (large traffic) and the Botanical Garden of Changchun University of Technology.
- They are analyzed by DP-LIBS: dual-pulse pre-ablation LIBS.
- A 1064 nm laser, 100 mJ, is focused on the surface of the soil sample. A pre-ablation 532 nm laser pulse is focused at a distance of 1 mm above the surface of the soil sample and some deviation of the focal point.
- The time delay between two pulses has influence on the enhancement of the spectrum in the process of DP-LIBS.




Double Pulse LIBS Experimental System.



Orthogonal pre-ablation DP LIBS configuration.




Influence of the pulse delay



Spectral enhancement.

- Device for identification and clasification of light materials in recycling processes.
- It combines a 3D detection system for the identification of the metallic pieces with the LIBS system for identification of their composition.
- It can support a flow of material of 4 tons/hour, meaning an analysis of ~100 particles per second with the LIBS system.
- A hybrid laser was used to perform this task: DPSS with flashlamp amplification.



Schematic structure consisting in a 3D scanner, the LIBS unit and compressed air nozzle.

- The integrated laser is responsible of two different tasks: first, there is a cleaning process of the surface; then, the spectrometric analysis is performed.


- CuZn-layer on Al 1200 alloy: 68  $\mu\text{m}$

A: Double pulse only

→ surface layer not penetrated


B: Tailored laser burst

→ bulk material is exposed by prepulses to double pulse LIBS analysis!



DPSS with flashlamp amplification:

1. Good beam profile, to focus in a small area of the sample material.
2. Double pulse generation, individual pulses well defined in shape and duration.
3. High repetition rate to allow the generation of the double pulse.
4. High energy in order to perform both tasks, cleaning and plasma generation for LIBS analysis.
5. The internal design of the laser must be as robust as possible to withstand the industrial environment (thermal and mechanical stresses, dust...).



- SpitLight 600 Laser with tailored burst.
- Detection and sorting rate 40 Hz
- Up to 600 mJ burst energy
- Moving samples at  $v= 3$  m/s
- Pilot demonstration for Al sorting  $> 4$  t/h
- In cooperation with




Foto: LIBS-sorting system at TITech GmbH, Mülheim-Kärlich, Germany

- Set-up installed at the DLR (German Aerospace Center), Institute of Technical Physics Langer Grund in Hardthausen (Germany)
- The DLR laser test range is designed for a wide field of laser application studies under central European atmospheric conditions: detection systems of biological, chemical and explosive hazardous substances.
- LIBS has been introduced for investigation of surface contaminants at distances up to 135 m.
- Micrometeorological measurements were performed simultaneously for a detailed description of the atmospheric properties affecting the propagation of the laser radiation.




Laser test range at DLR with LIBS set-up.

- Set-up installed at the DLR (German Aerospace Center), Institute of Technical Physics Langer Grund in Hardthausen (Germany)
- The DLR laser test range is designed for a wide field of laser application studies under central European atmospheric conditions: detection systems of biological, chemical and explosive hazardous substances.
- LIBS has been introduced for investigation of surface contaminants at distances up to 135 m.
- Micrometeorological measurements were performed simultaneously for a detailed description of the atmospheric properties affecting the propagation of the laser radiation.




Calculated beam spread due to optical turbulence after 135 m of free propagation.




Seasonal and diurnal influence on the optical turbulence.

- For plasma production a Nd:YAG laser is used: repetition rate of 10 Hz, maximum output energy of 800 mJ (@ 1064 nm) and pulse width of 8 ns.
- The first three harmonics of the Nd:YAG laser are used for creating the plasma.
- The laser beam is expanded by a Cassegrain telescope to a diameter of up to 290 mm and is focused to the target which is positioned at distances of 15, 35, 55, or 135 m, respectively.
- The light emitted from the generated plasma is collected by a 400 mm diameter Newton telescope and focused into an optical fiber with a diameter of 200  $\mu\text{m}$ .
- The collecting telescope covers an area with a diameter in the centimeter range.
- The detection is performed with a CCD spectrometer, triggered by the laser to ensure a proper timing (resolution of 0.3 nm and integration times from 1.28  $\mu\text{s}$  to 1.2 ms).




Schematically drawing of the LIBS set-up.

- To demonstrate LIBS capability for detection of explosives, black powder was used as a first representative: a mixture of sulfur, charcoal, and potassium nitrate (KNO<sub>3</sub>).
- Measurements performed at a bright sunny day.
- The characteristic doublet at 692 nm from KNO<sub>3</sub> was selected for black powder identification.
- Spectra at short distances  $\Rightarrow$  single pulse.  
Distances to more than 100 m  $\Rightarrow$  several laser pulses.
- More intense laser pulses compensate the drop of the detection signal with the increasing distance.
- Using adaptive optics increases also the quality of the measurements, as they compensate for atmospheric turbulence.




LIBS spectra of KNO<sub>3</sub> and black powder measured at a distance of 15 m.

- To demonstrate LIBS capability for detection of explosives, black powder was used as a first representative: a mixture of sulfur, charcoal, and potassium nitrate (KNO<sub>3</sub>).
- Measurements performed at a bright sunny day.
- The characteristic doublet at 692 nm from KNO<sub>3</sub> was selected for black powder identification.
- Spectra at short distances  $\Rightarrow$  single pulse.  
Distances to more than 100 m  $\Rightarrow$  several laser pulses.
- More intense laser pulses compensate the drop of the detection signal with the increasing distance.
- Using adaptive optics increases also the quality of the measurements, as they compensate for atmospheric turbulence.



- To demonstrate LIBS capability for detection of explosives, black powder was used as a first representative: a mixture of sulfur, charcoal, and potassium nitrate (KNO<sub>3</sub>).
- Measurements performed at a bright sunny day.
- The characteristic doublet at 692 nm from KNO<sub>3</sub> was selected for black powder identification.
- Spectra at short distances  $\Rightarrow$  single pulse.  
Distances to more than 100 m  $\Rightarrow$  several laser pulses.
- More intense laser pulses compensate the drop of the detection signal with the increasing distance.
- Using adaptive optics increases also the quality of the measurements, as they compensate for atmospheric turbulence.



Measured intensities of plasma emission (693.6 nm) at different distances (single pulse).



 GRUPO ALAVA

 Álava Ingenieros  
GRUPO ÁLAVA

 Preditec  
GRUPO ÁLAVA

 MRA  
GRUPO ÁLAVA

 Trazadia  
GRUPO ÁLAVA